

INFUSION SETS FOR UROLOGY

The infusion sets for urology have been specially developed to be as close as possible to the needs of urological washing.

One or two-way infusion sets are available, in particular the latter is allowing to use two bags consecutively, thus improving the efficiency of the operation.

The catheter cone fitting is universal and allowing connecting the infusion set to any type of catheter available


in the market; **infusion set with male Luer lock terminal fitting is also available.**

All infusion sets are equipped with drip chamber and roller to control and regulate the flow.

All infusion sets are equipped with drip chamber and roller to control and regulate the flow. Infusion units, UBH and Cystoflow bags create a synergistic mix of high quality to provide innovative solutions capable of giving a concrete response to the needs and requirements of the best urology departments.



INFUSION SETS FOR BLADDER WASHING

CODE	TYPE	TUBE		DRIp CHAMBER	CONNECTIONS	FLOW REGULATOR	
		LENGTh mm	Ø mm				
M060760S	1 VIA	1600	4,8x6,8	Ø 19 mm	STEPPED ENDPOINT	Roller	50
M060761S	1 VIA	2500	4,8x6,8	Ø 19 mm	STEPPED ENDPOINT	Roller	50
M060770S	2 VIE	1600	4,8x6,8	Ø 19 mm	STEPPED ENDPOINT	Roller	40
M060771S	2 VIE	2500	4,8x6,8	Ø 19 mm	STEPPED ENDPOINT	Roller	40
M060593S	2 VIE	2500	4,8x6,8	Ø 19 mm	L.L.M.	Roller	40

ALL CODES AND CHARACTERISTICS FROM PAGE 16