



# PUMP ADAPTERS

The use of administration sets addressed to specific infusion pumps is allowed by the addition of two- or


four-way pump adapter. They may be clear or amber coloured.



## TERMINAL CONNECTOR


Twist off to connection to the pump administrator spike.

### ADAPTERS IN PVC DEHP FREE

CODE	WITH N.R.V.	SPIKE	CONNECTIONS	
MCI/350A	2	2 WAY- AIR FILTER 0,2 µm	SPIKE MEMBRANE	50
MCI/350B	4	2 WAY- AIR FILTER 0,2 µm	SPIKE MEMBRANE	50

## PHOTOSENSITIVE DRUGS

### ADAPTERS IN PVC DEHP FREE FOR PHOTSENSITIVE DRUGS

CODE	WITH N.R.V.	SPIKE	CONNECTIONS	
MCI/350AA	2	2 WAY- AIR FILTER 1,2 µm	SPIKE MEMBRANE	50
MCI/350BA	4	2 WAY- AIR FILTER 1,2 µm	SPIKE MEMBRANE	50