

DIMAR

PRODUCT LINE

Intensive Care Unit &
Hospital Departments

100%
MADE IN ITALY

DimAir Line

FIFTH
FIXATION POINT

INTERCHANGEABLE
HEADGEAR

VERSIONS
-SINGLE PATIENT
- RESTERILIZABLE

MORE SPACE FOR
THE NOSE

FIBROBRONCHOSCOPE
DEDICATED PASSAGE

INTERCHANGEABLE
CONNECTORS

NASOGASTRIC
TUBE
DEDICATED
PASSAGE

TOTAL FACE DIMAX

DIMAX TOTAL FACE MASK

Total Face mask for non-invasive mechanical ventilation in PSV, BIPAP and CPAP techniques. Featuring 5 attachment points, it optimizes patient/ventilator synchronization and pressure distribution to the face, reduces leakages and improves comfort. Increased space for prominent noses, combined with the ability to change connectors and switch the same mask between vented and non-vented modalities, makes DIMAX the ideal interface for non-invasive ventilation.

7 plus

- + Connection
- + Comfort
- + Sealing
- + Practicality
- + Flexibility
- + Safety
- + Versions



The 5th fastening point allows greater grip on the patient's forehead



Kit with Easy Vent (700/3482 Filtered PEEP)



"Non Vented" connector with anti asphyxia valve

TECHNICAL FEATURES:

- 5 fixation points.
- More space for the nose.
- Dedicated access for Nasogastric tubes
- Interchangeable Connectors:
 - 22M Vented
 - 22F Non Vented



Dedicated passage for NGT



"Non Vented" Connector for bitube ventilator



"Vented" connector with anti asphyxia valve for single tube ventilator

ADVANTAGES:

- Better patient/ventilator synchronization.
- Improved pneumatic seal.
- More comfort for prominent noses.
- Lower headgear traction.
- Reduced risk of contact injury.
- Possibility of passage for NGT and probe for Fibrobronchoscopy, without leakage and interruption of ventilation
- Single-use and resterilisable versions in autoclave 134° C. - Gas Plasma - ETO



REFERENCES:

SINGLE PATIENT *	SIZE	RDM	CND/EMDN
700/6989	XS/PED	2683312	R0301010202
700/6709	SMALL		
700/6710	MEDIUM		
700/6711	LARGE		
700/6712	EXTRA LARGE		
700/7243	EXTRA EXTRA LARGE		

RESTERILIZABLE*	SIZE	RDM	CND/EMDN
700/7268	XS/PED	2683333	R0301010202
700/6739	SMALL		
700/6740	MEDIUM		
700/6741	LARGE		
700/6742	EXTRA LARGE		
700/7269	EXTRA EXTRA LARGE		

* In single pack, multipack of 5 pieces

** In single pack



HMC Premedical S.p.A

Via Bosco 1/3 – 41037 Mirandola – MO – I
Tel: +39 0535 22704 Fax: +39 0535 609546

