

DIMAR

PRODUCT LINE

Intensive Care Unit

DimAir Line



CPAP HELMETS
PEDIATRIC AND NEONATAL

PEDIATRIC AND NEONATAL CPAP HELMETS

Single-patient helmets for the treatment of respiratory failure in CPAP and BI-PAP techniques in neonate patients from 3 to 5 kg and pediatric patients from 5 to 15 kg. Equipped with a zip closure, they allow quick and easy access to the patient. Convenient and practical.

TECHNICAL FEATURES:

- Airtight closure with zipper.
- Made of biocompatible plastic material.
- Symmetrical Right / Left connections for convenient and easy connection to the equipment on either side.
- Hygienic antibacterial containment flap Coolmax®.
- Sealed access from 3 to 19 mm.
- 2 fastening eyelets on the top of the helmet for an alternative/ additional fixing.
- Connection ISO 22M, 22F.
- Latex and PHT free



The "eyelet" system allows additional fastening.

ADVANTAGES:

- Quick access to the patient's face.
- Possibility to perform nursing, feeding and weaning without removing the helmet.
- Extremely comfortable, practical and effective.
- Avoids oro-nasal tracheal intubation infections.
- Without the rigid ring, it avoids decubitus to the patient's nape.
- The gas inlet and outlet can be reversed or connected on the same side; so it is also possible to place the patient on one side.
- Prevents the formation of necrosis or pressure sores caused by facial masks or nasal tubes.
- In combination with the DIMAR PEEP valve, it ensures bidirectional flow in case of lack of pressure, preventing the risk of asphyxia.

REFERENCES:

CODE	SIZE	RDM	CND/EMDN
500/5447	NEO	2683503	R0301010401
500/5448	PED		

- In single pack
- Multiple pack of 5