

DIMAR

PRODUCT
LINE

Ambulance
Emergency
Hospital Department

DimAir Line



**COMFORT COULISSE
CPAP HELMET**


EASY VENT CPAP HELMET EMERGENCY COULISSE

EMERGENCY COULISSE CPAP HELMET (AMBULANCE, EMERGENCY, FIRST-AID)

- Disposable helmet for emergency treatment of acute respiratory failure in CPAP. Equipped with quick closure "Drawstring", can be quickly applied by a single operator.
- One size fits all, ensures airtight seal on patients of any size.
- Provided in small packages for reduced ambulance footprint.
- Equipped with high efficiency Venturi Easy Vent generator, adjustable PEEP valve and pressure control manometer, it can be connected directly to the oxygen cylinder for the administration of CPAP in patients with cardiogenic edema directly from the patient's home.
- Excellent quality / price ratio.

TECHNICAL FEATURES:

- One size fits all.
- Drawstring closure.
- Venturi generator, PEEP valve and pressure gauge pressure already inserted for quick use.
- Soft and flexible armpit straps for quick locking.
- Made of biocompatible plastic.
- Single patient.
- Latex and PHT free.



ADVANTAGES:

- Quick positioning, opening and closing.
- Can be positioned by a single operator.
- Rapid pressurization of the respiratory tract.
- Applicable to patients of any size.
- Compact "pocket" package for reduced footprint, can be placed in any space.
- Practical and effective.

REFERENCES:

CODE	SIZE	Internal volume* (Lt – Head excluded)	RDM	CND/EMDN
500/8612I	ONE SIZE FIT ALL	15,7 Lt	2683503	R0301010401

- Single bag
- Multiple pack of 5
- Possibility to create custom models, versions and packages on request

*+/- dimensional tolerances excluding collar and patient head



HMC Premedical S.p.A

Via Bosco 1/3 – 41037 Mirandola – MO – I
Tel: +39 0535 22704 Fax: +39 0535 609546

