

DIMAR

PRODUCT LINE

Ambulance
Emergency
& Hospital Departments

DimAir Line



***HIGH PERFORMANCE
MAX EASY VENT***

High performance venturi system, for CPAP therapy.

Allows selection of FiO₂ values, predefined at 30% - 60% - 90% .
Pre-assembled and ready-to-use device, with reduced weight and minimal footprint, it offers easy handling, practicality and speed of execution for the operator and effectiveness for the patient.

HIGH PERFORMANCE MAX EASY VENT is the practical and effective solution for administering air/oxygen mixtures in CPAP; Easy and reliable compact device in Emergency, Urgency and Therapy.

TECHNICAL FEATURES:

- Corrugated tube 180 cm
- UNI, DIN, AFNOR quick couplings
- Pressure monitoring manometer



Total Face Mask



ADVANTAGES

- Adjustable
- Low weight, compact
- Ready for use



Oronasal Mask



DESCRIPTION	CODE	SIZE	PACKAGING	RDM	CND/EMDN
HIGH PERFORMANCE MAX EASY VENT <i>With Oronasal Mask</i>	700/7580	S	In single pack. Multipack of 10 pieces.	R0301010202	2683312
	700/7581	M			
	700/7582	L/XL			
HIGH PERFORMANCE MAX EASY VENT <i>With DIMAX Mask</i>	700/7621	S	In single pack. Multipack of 5 pieces.		
	700/7622	M			
	700/7623	L			
	700/7624	XL			
	700/7625	XXL			

